
#### NEW EQUIPMENT

2018 Mid-Year Progress Report



### NEW EQUIPMENT

### 5.11TACTICAL BLS RESPONDER BAG

 Replaces BLS bag on Tanker and Brush Truck

Meets County standards

Standardized layout for both bags



FIRE

ALS

BLS

# MASS CASUALTY RESPONSE GEAR

 Rapid Triage and Life-Saving Treatment Kit



FIRE ALS BLS

## MASS CASUALTY RESPONSE GEAR

- Still in development (target roll-out date August 2018).
- Medic 79 and LS 7 will carry 2
- Engine 72 will carry I
- An additional "First 15-Minute Kit" will be developed and deployed on Medic 79 at a later date.
- Specific and required training program will be provided prior to deployment.

#### IPAD & KEYBOARD

- Apple iPad 9.7"
- ZAGG Bluetooth Keyboard & Case



FIRE ALS BLS

#### IPAD & KEYBOARD

- Though ruggedized, neither the case nor iPad are water resistant. DO NOT EXPOSE TO MOISTURE OR WATER.
- To be used for completion of FireManager Forms and EMEDS Elite ePCRs (After 9/1/18).
- Replaces the Toughbook with 9.7" iPad
- Rugged case with detachable, bluetooth keyboard.
- · Backlit keyboard can be used in low-light situations

# LUCAS 3 BATTERY CHARGER

 Desktop battery charger for LUCAS 3 batteries



ALS


BLS

### LUCAS 3 BATTERY CHARGER

- Will allow for the LUCAS to be stored more appropriately.
- Batteries should be swapped during the weekly unit check.
- Use caution removing batteries, lifting up slightly before removing from charger, to prevent severe damage to battery.
- Battery charger to be stored in "stack" cabinet, below MedVault.

#### ORALTHERMOMETERS

- Hospital grade oral thermometer
- Stored in the Action Area cabinet on Medic 79
- Stored in the main section of the BLS Bag on LS 7.







# T-POD PELVIC BINDING DEVICE

 Disposable device for stabilizing pelvic fractures.





# T-POD PELVIC BINDING DEVICE

- · Trauma Pelvic Orthotic Device for pelvic fracture stabilization.
- Will be carried in the Severe Trauma Compartment in M 79.
- Will be carried in the new Frac-Pack on LS 7.
- All providers are required to complete the online training program at www.lvfdems.com.
- It is no longer acceptable to use the KED Board for pelvic binding.

# ACR-4 AMBULANCE CHILD RESTRAINT SYSTEM

 Crash-test pediatric transport device





# ACR-4 AMBULANCE CHILD RESTRAINT SYSTEM

- · Safely secure pediatric patients of all sizes using one system.
- System allows for full access to patient for enhance safety and optimal patient care.
- Specific storage location on Medic 79 still TBD.
- All providers are required to complete the online training program at www.lvfdems.com.

#### PRE-LOADED ETTUBES

- Cuffed Endotracheal Tubes will now come with a stylet preloaded.
- ETT Sizes:
  - 6.0, 7.0, 8.0, 9.0
- Not to be used with Video Laryngoscope





### AIRTRAQ A-390

Video Laryngoscope System





### AIRTRAQ A-390

- Can be used on Pediatric and Adult intubations
- Smaller profile, bigger screen
- Failsafe feature integrated into blade

Size / Color	ET Tube Size
0 - GREY	2.5 - 3.5
I - PURPLE	4.0 - 5.5
2 - GREEN	6.0 - 7.5
3 - BLUE	7.0 - 8.5

- Integrated video recording, no wires or cables required
- Will be carried on Medic 79 and LS 7, specific location to be determined
- All ALS providers are required to complete the online training program at www.lvfdems.com

#### KING LTS-D TUBES

 Rescue airway for ALS Providers





#### KING LTS-DTUBES

- King LTS-D Tubes are approved for patient's of all sizes.
- · I complete set will be stored in the Airway cabinet in Medic 79.
- Sizes 0, 1, 2, and 2.5 will be stored in the intubation pack of each Pediatric Bag
- Sizes 3, 4, and 5 will remain in the front pocket of the red ALS Bag

Size / Color	Weight
0 - TRANSLUCENT	< 5kg
I - WHITE	5 - 12 kg
2 - GREEN	12 - 25 kg
2.5 - ORANGE	25 - 35 kg
3 - YELLOW	4 - 5 Ft
4 - RED	5 - 6 Ft
5 - PURPLE	> 6 Ft

### QUESTIONS?

See an EMS Officer

Matthew Burgan, EMS 7 mburgan@LineboroVFD.org 410-365-2610 Summer 2018